

**12 O'CLOCK LEFT/RIGHT WORKSTATION
WITH 3 HP AUTOMATIC DELIVERY & ASSISTANT'S PACKAGE
PLUMBING DIAGRAM**



Tubing Legend			4 Hole Handpiece Tubing Legend		2 Hole Syringe/Signal Tubing Legend	
Part#	Reference	Description	Reference	Description	Reference	Description
1203		Supply Tubing, 1/8", Poly Clear		Drive Air		Ribbed
1204		Supply Tubing, 1/8", Poly Green		Air Coolant		Non-Ribbed
1205		Supply Tubing, 1/8", Poly Gray		Water Coolant		
1206		Supply Tubing, 1/8", Poly Orange	N/A	Exhaust		
1207		Supply Tubing, 1/8", Poly Red				
1208		Supply Tubing, 1/8", Poly Yellow				
1402		Supply Tubing, 1/4", Poly Blue				
1405		Supply Tubing, 1/4", Poly Gray				
1407		Supply Tubing, 1/4", Poly Red				

