

*DCI* edge



DCI EDGE SERIES 5 OPERATORY EQUIPMENT

# OUR VALUES, OUR STORY

## Crafted with Decades of Experience

John and Janelle Spencer had a dream to start their own business. They knew that they had the discipline, drive, and passion to create a successful business and give back to their community. What began with a small home loan, an idea for a syringe, and a makeshift office in their living room grew into a global parts giant.

With the help of their children, Jason and Amy, the Spencers continue to change the face of dental with continuous innovation, remarkable expansion and a dedication to the highest standards. To this day, they remain true to their core values, making their customers' lives easier.

Jason, Amy & John Spencer



*Newberg, Oregon*

## The Pacific Northwest

DCI first opened its doors in Newberg, Oregon, the heart of dental equipment manufacturing in America. Over 350 employees are part of the DCI family which has proudly served the dental industry over 35 years.



*Charlotte, North Carolina*

## East Coast Expansion

In 2018, DCI opened an additional 40,000 sq ft cabinet manufacturing facility. Equipped with state of the art CNC machinery, DCI manufactures cabinets on both coasts.



## Reliable Inside & Out

DCI is the largest manufacturer of dental parts in the world. We supply most equipment manufacturers in the market with delivery system parts – the parts that actually make your unit work, the “engine” of your dental system.

Most likely, you’ve been using DCI parts for years without even knowing it. From foot controls to control blocks, vacuum valves to quick disconnects, tubing to flow control valves, there’s a good chance DCI is inside the units you’ve grown accustomed to over the years.




1983	1991	1993	1997-2005	2016	2018
<b>Company Beginnings</b> DCI is founded by John & Janelle Spencer.	<b>New Delivery System</b> DCI designs its first world acclaimed dental system.	<b>Global Growth</b> DCI becomes largest dental parts manufacturer in the world.	<b>Expansion</b> DCI expands through the acquisition of Pelton & Crane and Marus, building them into the second largest dental company in the US. Companies are sold in 2005.	<b>Introducing DCI Edge</b> DCI returns to the Equipment market with the introduction of the new Series 5 Operator Package.	<b>New Product Innovation</b> DCI continues to expand its dental offering with the introduction of a full cabinetry and orthodontic line.

# SMART INVESTMENT

## Why Pay More?

DCI Edge products offer the features you want, the reliability you need, and the comfort your patients deserve for significantly less money than other popular brands.



*DCI Edge*

## Financial Freedom

It's your practice. It's your budget. DCI Edge delivers quality equipment at a price point that allows you to truly enjoy life.



## Reinvest

Upgrade your practice with the latest & greatest technology and secure your edge in the ever changing tech world.



## Save

There is a peace of mind that comes with money in the bank. Save for a rainy day fund or your retirement plan.




## Enjoy

Buy your dream car. Play the course you have always wanted. Experience a new destination from your bucket list.

# BEAUTIFUL INSIDE AND OUT



**Built to Last**  
Our design philosophy blends innovation with high quality construction. We are proud to show you what's under the hood.



Cutting Edge diagnostic capabilities that displays data from your chair on a personalized App.

# GET COMFORTABLE STAY HEALTHY

## Access to the Oral Cavity

Unique contoured back design ensures full access between the 3 and 9 o'clock working positions while comfortably supporting the patient.




## Ergonomics

Thin, narrow back and dual articulating headrest allows for ideal positioning for the provider and patient. This leads to enhanced comfort and efficiency for the entire dental team.




## How Low Can You Go?

The Series 5 chair height ranges from 16" to 31" at the seat so you can comfortably practice sitting or standing.



## Are you in the right position?

Optimal positioning near the patient ensures a healthy posture, reducing pain and fatigue.




# KEEP YOUR PATIENTS SMILING



DCI *EDGE*

## Multi Layer System

Four different layers of foam with varying levels of density deliver premium comfort and support. Adapts to all body types.



## Soft Start and Stop

Adjust your chair without disrupting your patient with smooth start and stop technology.



## Patient Support

The innovative headrest scoop and back design provides patients with back support regardless of height.

## Design Appeal

Available in wide or narrow back, plush or asepsis, in a large variety of colors. Upholstered armrests add to the comfort and style.

## Doctor & Assistant Stool

Advanced seat ergonomics: tilt and height adjustment add to overall comfort.



ErgoBack Plush

NarrowBack Asepsis

NarrowBack Plush

ErgoBack Asepsis



# PERFORMANCE MEETS DESIGN

## Ergonomic Access to Instruments

Additional head pivot swivels to improve access & reduce unnecessary movement.



## Internal Components

Reliable pneumatic components & ample room for high tech integration.

## Dual Integrated Touchpads & Brake Handles



Lights On/Off



Return To Last Position



Programmable Buttons for Preset Positions



Backrest Incline and Recline



Chair Up and Down



## Advanced Delivery System

DCI's Series 5 delivery system integrates technology to improve operator efficiency and eliminate workspace clutter. Enjoy the flexibility of integrating the devices you use today as well as what you might need in the future.

# DENTAL SYSTEM SOLUTIONS

## Swing Mount Delivery System

Delivery system and light easily convert from left to right handed operation. The compact swing radius optimizes space making it perfect for operatories of any size.



Shown with Rear Cuspidor

## Post Mount Delivery System

The post mount delivery system features an integrated vacuum trap, water system, as well as air and water quick connects.



Shown with Cuspidor



# REFLECTIVE LED

## Low Cost of Ownership



### Energy Efficiency

Requires 80% less power compared to traditional dental lights.



### Longer Bulb Life

Significant Savings on Replacement Bulbs: 10-15 year expected life vs. annual replacement.

## Patient and Operator Comfort




### Cool Running

LED technology eliminates heat output, maximizing patient and operator comfort.



### Precise Illumination

Prevents patient glare.



*DCI edge*

## Less Eye Fatigue Increased Visibility

Optimal light performance reduces shadows, improving illumination of the oral cavity.


## Uniform Light Pattern

DCI delivers cutting edge illumination technology by utilizing over 500 reflective facets to bounce light back to the oral cavity. The result is a uniform light pattern that increases visibility for doctors and dental assistants.

# BRILLIANT ILLUMINATION



## Integrated Touchpads



## Mounting Configurations



## No-Touch Sensors

Control the light without contact using programmable no-touch sensors.



## No-Cure Setting

Filters out UV light to prolong the working time of photo-initiated resins.



# Dental Cabinetry Should Be

**Functional.  
Beautiful.  
Durable.  
Streamlined.  
Affordable.**





Designed  
with **your**  
**needs** in  
mind

### Get Comfortable

Ergonomic operatory design allows for simple class I, II, III movements, eliminating the need to twist or turn.

### Save Time

Multi-task without procedure interruption using pass-through doors. Front and side configurations keep items easily accessible.

### Stay Connected

Cleanly connect all of your ancillary products using pre-wired USB and GFI outlets. Fully functional clock timer enhances productivity.

### Expanded Workspace

Maximize countertop space so you have adequate room to practice efficiently. Center console hides cords and cables.



FUNCTIONAL

# BEAUTIFUL

## Signature Design Aesthetic

Color options curated by cutting edge designers.

## Focused on the Details

Patented beveled Infinity Edge countertops add an elegant look to your practice.

## A Perfect Finish

1 mm Fusion edgebanding elevates the finish of your cabinet.

## Quiet Environment

Soft close doors and drawers create a peaceful patient environment.

Take  
**pride** in  
your practice



**Built**  
with care  
& precision.



### High Quality Construction

Heavy duty structural steel support and base are built to last.

### Precision manufacturing

Designed using industry leading software & cut using state of the art CNC Machinery, ensuring quality at the highest level.

### Reliable Components

Engineered with top of the line components to stand the test of time.

### Designed for Dental

Constructed to support all dental operator needs including units and lights.




# STREAMLINED




## A Workspace Built Around You

Sterilization is one of the most important steps in a practice, and one size never fits all. With DCI Edge, you can create a Sterilization Center that fits the needs of your practice. All options feature a process based flow, to keep your productivity high and your risk extremely low.


Straight Run




Galley Style



L Shape



U Shape



Receiving, Cleaning & Decontamination



Preparation & Packaging



Sterilization



Sterile Storage & Make Ready

### Sterilize with Standards

Designed to follow CDC guidelines with exceptional focus on infection control and patient safety.



### Sterilize with Style

Fully functional clock timer enhances productivity, while stainless steel valance adds to the beautiful design aesthetic.



# Sterilization SIMPLIFIED



## Flip Down Wrap Station

Two conveniently spaced drawers designed to pull out & flip down ensuring ample space to wrap any size cassette.



*DCI* **EDGE**

## Systematic Decontamination

Increase flow with specialized waste containers, stainless steel undermount sink, dual integrated drainboards & conveniently placed glove dispensers.



## Factory Integrated Accessories

Our designers and engineers have curated a list of the top options in dental sterilization. Factory installation provides optimal placement and smooth workflow.



## Sterilization Towers That Are Built To Last

Built with water proof materials to guarantee protection from moisture damage. Multiple tower options accommodate your specific sterilizer selections.



## Simple Sterile Organization

Pull out storage with wire racking can hold numerous cassettes, tubs, trays, & wrapped instruments. Multiple adjustable shelves allow for fast and easy management.



## Clear Visual Workflow

Red Lights indicate contaminated storage, while blue lights illuminate sterile storage. Additional task lighting brightens the entire worksurface for maximum visibility.



# Find the Perfect Color

## Over 85 designer choices


Upholstery



Premium Laminates



Solid Surface



For our full color selection, visit us at  
[www.dciedge.com/request-samples](http://www.dciedge.com/request-samples)

