

12 O'CLOCK FLOOR-MOUNTED DELIVERY SYSTEMS



The DCI Edge Floor-Mounted Delivery Unit enables ergonomic access to instruments, offering vertical adjustment from left to right. The floor-mounted design is robust and intelligently crafted to neatly fit into the cabinet, delivering strong support, heightened durability, efficient space utilization, ample clearance for stool casters, and an enhanced aesthetic appeal.



DELIVERY UNIT FEATURES

- Improved floor-mount design provides 6" side-to-side movement without encroaching on operator space.
- Accommodates ancillary integration, now or in the future.
- Mobile work surface has an adjustable height range of 3".
- Unique design accommodates easy positioning for right- or left-handed practitioners.
- Alignment of collection trap and water bottle enable convenient service position for Assistant.
- Assistant touchpad pivots 180 degrees.
- Work surface tray holder provides additional space for supplies.



12 O'CLOCK FLOOR-MOUNTED DELIVERY SYSTEMS

	Part Number		Retail Price
12 o'clock Floor-Mounted Delivery System with Dr. & Assistant's Instruments, Mobile Solid Surface Work Top & Touch Pads Can only be paired w/DCI Edge Series Chair	DR5550	Gray	\$11,000.00
12 o'clock Floor-Mounted Delivery System with Dr. & Assistant's Instruments, Mobile Solid Surface Work Top	DR5450	Gray	\$10,500.00
12 o'clock Floor-Mounted Assistant's Instruments, Mobile Solid Surface Work Top & Touch Pads Can only be paired w/DCI Edge Series Chair	DR5551	Gray	\$4,950.00
12 o'clock Floor-Mounted Assistant's Instruments & Mobile Solid Surface Work Top	DR5451	Gray	\$4,750.00
12 o'clock Floor-Mounted Solid Surface Work Top Only	DR5454	Gray	\$3,500.00



42" 12 O'CLOCK CABINET with 12 o'clock Floor-Mounted Dr. & Assistant's Delivery System



DR5551 Floor-Mounted Assistant's Delivery

Floor-Mounted Delivery standard features include:

- Solid Surface Work Top
- Mobile Work Surface Tray
- Air and Water Quick Disconnects

