

DCI edge



Our values, **OUR STORY**

PORTLAND OREGON, 1982

John and Janelle Spencer had a dream to start their own business. They knew that they had the discipline, drive, and passion to create a successful company and give back to their community. What began with a small home loan, an idea for a syringe, and a makeshift office in their living room grew into a global parts giant.

1983

COMPANY BEGINNINGS

DCI is founded by John & Janelle Spencer.

1991

NEW DELIVERY SYSTEM

DCI designs its first world acclaimed dental system.

1993

GLOBAL GROWTH

DCI becomes largest dental parts manufacturer in the world.

1997-
2005

EXPANSION

DCI expands through the acquisition of Pelton & Crane and Marus, building them into the second largest dental company in the US. Companies are sold in 2005.

2016

INTRODUCINGDCIEDGE

DCI returns to the Equipment market with the introduction of the new Series 5 Operator Package.

2018

NEW PRODUCT INNOVATION

DCI continues to expand its dental offering with the introduction of a full cabinetry and orthodontic line.

CRAFTED WITH DECADES OF EXPERIENCE

With the help of their children, Jason and Amy, the Spencer's continue to change the face of dental with continuous innovation, remarkable expansion and a dedication to the highest standards. To this day, they remain true to their core value, making their customers lives easier.

"Sometimes all it takes is an idea and passion to change the course of your life and career."

DCI *EDGE*



West Coast Headquarters

DCI first opened its doors in Newberg, Oregon, the heart of dental equipment manufacturing in America. Over 350 employees are part of the DCI family which has proudly served the dental industry for over 37 years.



East Coast Headquarters

In 2018, DCI opened an additional 40,000 sq ft cabinet manufacturing facility in Charlotte, North Carolina. Equipped with state of the art CNC machinery, DCI manufactures cabinets on both coasts.



Ortho CHAIRS



Fixed Base

Equipped with a solid steel base to ensure maximum stability. The Integrated foot control makes positioning your patients quick and easy.



Adjustable Base

Height ranges from 16" to 31" at the seat so you can comfortably practice sitting or standing. Available with multiple backrest and upholstery options.



V HEADREST



PONYTAIL HEADREST



DUAL ARTICULATING
HEADREST

DCI **EDGE**



Fixed & Adjustable Chairs
are available with optional
armrests.



Ortho CABINETS



RECESSED GROOVE

Prevents orthodontic bands and brackets from rolling off of the countertop.



ADVANCED FEATURES

Soft close doors and drawers. Multiple configurations available. Premium high pressure laminate and solid surface countertop options.

HEAVY-DUTY CONSTRUCTION

High quality steel and aluminum castings stand the test of time. Fully hydraulic lift and back mechanisms for high powered weight support.

Standard Mobile



Narrow Fixed Base



Narrow Mobile



Standard Fixed Base

Seamless **INTEGRATION**

Delivery Module

High quality components
for precision and speed.



CPU Package

Computer module-flat screen monitor mount includes CPU cooling fan.



*Monitor not included



Quick Disconnect Panel

Simple & clean installation when adding accessories.



Glove Dispenser

Integrated glove dispenser to reduce clutter.



Integrated Trash Receptacle

Easy access for quick and efficient clean up.

Reflective LED TECHNOLOGY

NO-TOUCH SENSORS

Control the light without contact using programmable, no-touch sensors.



DURABLE DESIGN

DCI's heavy-duty metal structure and handle design along with high quality components demonstrate a strong commitment to building the highest quality dental equipment.



NO-CURE SETTING

Filters out UV light to prolong the working time of photo-initiated resins.

Ceiling Mount Light



Track Mount Light



Wall Mount Light



Post Mount Light



LESS EYE FATIGUE INCREASED VISIBILITY

Optimal light performance reduces shadows, improving illumination of the oral cavity.

UNIFORM LIGHT PATTERN

DCI delivers cutting edge illumination technology by utilizing over 500 reflective facets to bounce light back to the oral cavity. The result is a uniform light pattern that increases visibility for doctors and dental assistants.

ACCURATE SHADE MATCHING

5000 Kelvin color temperature with a color rendering index (CRI) of 90 produces a neutral white light, which is ideal for color matching at any intensity level.

FIND THE PERFECT COLOR

Over 85 designer choices


UPHOLSTERY



PREMIUM LAMINATES



SOLID SURFACE



For our full color selection and to order samples visit us at www.dcedge.com/request-samples

