



SERIES 5 DENTAL CHAIR
 Technical Specification

DCI edge



CHAIR WEIGHT: 300lbs / 136.4kg
 UPHOLSTERY WEIGHT: 30lbs / 13.6kg
 SHIPPING WEIGHT:

POWER REQUIREMENTS: 115VAC, 8A, 50-60Hz
 230VAC, 2.5A, 50-60Hz
 POWER CORD LENGTH: 5 Feet (1.5m)

IEC MEDICAL CLASSIFICATION I
 TYPE: B Chair upholstery
 OPERATION: Intermittent
 DUTY CYCLE: 20 Seconds ON, 5 Minutes OFF
 SPLASH PROTECTION: IPX0 (FOOT SWITCH IPX1)

CERTIFICATIONS: IEC 60601-1 3rd Edition,
 CAN/CSA C22.2 No. 601.1

