

DCI Edge

Orthodontic Equipment Solutions:

Performance + Value

Explore orthodontic equipment solutions designed for today's dental teams. Customize your setup for greater efficiency and functionality with ergonomic features that maximize your office space. Backed by a 10-year warranty, DCI Edge ortho equipment is a smart, reliable investment for your practice.



DCI Edge Ortho Chairs

Engineered for reliability and durability, our chairs are designed by experts who understand how you work. We've eliminated unnecessary, costly features to provide exactly what you need to run your practice efficiently—at a budget-friendly price.

Fixed Base

DCI Edge Orthodontic Chair

Features a solid steel, fixed base for maximum stability.

- Adjustable backrest that helps you position patients quickly and comfortably
- Integrated footswitch with preprogrammed positions
- Auto-exit makes patient positioning even faster and easier



Adjustable Base

DCI Edge Series 4 Chair

Offers a greater range of movement, allowing for both sitting and standing procedures.

- Seat height adjusts from 16" to 30"
- Dual articulating headrest
- Multiple backrest options



FLAT HEADREST



DUAL ARTICULATING HEADREST



V HEADREST



PONYTAIL HEADREST

Customize each chair with your choice of headrest style.

DCI Edge Ortho Cabinetry

Optimize your operatory with cabinetry designed for durability and flow. Featuring quiet, soft-close doors and customizable configurations, our high-pressure laminate cabinets are built to last. Includes solid-surface countertop options. A recessed countertop groove helps prevent orthodontic bands and brackets from rolling off, keeping your workspace organized and efficient.

Heavy-Duty Construction & Support

- Durable steel and aluminum castings built to stand the test of time
- Fully hydraulic lift and back mechanisms for powerful weight support
- Available with fixed or mobile base options
- Standard and narrow versions available
- Customizable drawer and door configurations



Free Standing Ortho Cabinet, Available in standard 28" and narrow 24"

Our fixed cart features a 6-inch steel subbase that protects your cabinetry from water damage and can be installed over existing utility hookups, keeping your space tidy and uncluttered.



Mobile Ortho Cabinet, Available in standard 28" and narrow 24"

Our mobile cart allows you to move your utility hookups with you, providing a more flexible option and allowing you to maximize your space with additional operatories.



Customize the Control Head Delivery Unit

High-quality components designed for precision and speed. Standard configuration includes HVE with tubing, SE with tubing, DCI syringe with tubing, and HP tubing (4-hole asepsis). Control head slides out easily when unlocked for quick, convenient maintenance.

Additional options include

- CPU upgrade
- Keyboard tray
- Mouse tray
- Dual glove dispenser
- Sanitizer wipe holder
- Waste chute/cutout

Boost care and comfort for young patients

Our patented design Kid Cushion is available in all upholstery colors—perfect for adding a fun pop of color or creating a seamless, coordinated look throughout your dental office.

