

DCI *Edge*

Series 4 Equipment Solutions

Blending Value and Style

DCI Edge Series 4 offers the perfect blend of value and style. With all the essential features your practice needs to run smoothly—at a price that makes sense for your budget.



DCI Edge Series 4 Chair

Our classic chair includes everything you need—at a budget-friendly price.

Patient and Practitioner Comfort

Thin, narrow back design provides comfortable, ergonomic access to the oral cavity while comfortably supporting the patient.

Easy Patient Access

Foldout armrests provide easy entry and exit.

Quick Headrest Positioning

Dual articulating headrest facilitates quick repositioning.

Heavy-Duty Construction

High-quality steel and aluminum castings stand the test of time. Fully hydraulic lift and back mechanisms for high-powered weight support.

Customized Comfort

Customize each chair with your choice of chair back and headrest styles.



FLAT HEADREST



DUAL ARTICULATING
HEADREST



V HEADREST



PONYTAIL HEADREST



ERGOBACK PLUSH



ERGOBACK ASEPSIS



NARROW BACK PLUSH



NARROW BACK ASEPSIS

DCI Edge Series 4 Delivery Systems

The Series 4 is the only delivery system that can retrofit most dental chairs. Configurable with up to five handpieces/instruments and one syringe. Each handpiece has individual drive air and water coolant adjustments.



Series 4 Swing Mount Delivery System

Delivery system and light easily convert from left- to right-handed operation for true ergonomic positioning.

Series 4 Hygiene Delivery

Integrated vacuum lines provide a sleek, modern look. Available in multiple configurations.

Series 4 Side Delivery

Available in either white or gray to coordinate with existing operatory equipment.

Series 4 Cart and Operatory Support Cart

Flexible solutions for even the most challenging operatories.



Series 4 Post Mount Delivery System

Features an integrated vacuum trap, water system, optional cuspidor, as well as air and water quick connect.

Additional Series 4 Equipment Options



LED LIGHT

Get reflective LED technology at an affordable price.



TELESCOPING ASSISTANT'S PACKAGE

Featuring ergonomic syringes and vacuum valves.



REAR ASSISTANT'S PACKAGE

Mounts to rear of chair for flexible positioning.



CUSPIDOR

Durable metal construction with porcelain bowl.

Ready to Level Up?

DCI Edge Series 5 provides a more premium patient and provider experience with a range of additional features.



Compare
DCI Edge Chairs